



Oil Soluble Copper Tripeptide-1 (Oil Soluble GHK-Cu)

Product information

Product: Oil Soluble Copper Tripeptide-1

INCI Name, Case# & Ratio

Copper Tripeptide-1 49557-75-7 0.1%

Glycerin 56-81-5 5.0%

1,2-Hexanediol 6920-22-5 2.0%

Triheptanoin 620-67-7 40.0%

C13-15 Alkane 64742-46-7 35.0%

Polyglyceryl-2 Isostearate 73296-86-3 13.0%

Aqua 7732-18-5 4.9%

Properties: Light blue clear liquid

PH: 3.6-6.0

Solubility: Good oil solubility

Specifications: 1000ppm, 2000ppm or custom ed

Package: 1.0Kgs/Bottle or 5.0kgs/Bottle



Patented Technology Highlights

On the basis of the advantages of small molecular weight, easy absorption and wide compatibility of peptide raw materials, the oil-soluble peptide solution strengthens the stability of the peptide and expands the application field to avoid the phenomenon that peptides are very sensitive to temperature and easy to change color, turbidity or delamination; Oil-soluble peptide solution, using ingredients with a natural origin index of 1, is characterized by oil-soluble active ingredients of

TECHNICAL DATA SHEETS

PROFESSIONAL PEPTIDE MANUFACTUR

emulsifiers and solvents, which is light and non-greasy, skin-friendly, and can be quickly absorbed by the skin, and has strong compatibility with common oils.



Functions & Usage

Promote the growth of fibroblasts and keratinocytes, protect injured tissues from further damage by oxidation, and promote the healing of damaged skin. It has the action of peptide GHK, which promotes the synthesis of collagen and elastin, tightens loose skin, and reduces wrinkles and deep wrinkles. It can be added to various kang wrinkle repair products and hair growth products such as face, eyes, neck, etc., such as functional liquid, nutritional cream, essence essence, mask and sunburn cream, BB cream.

Recommended Concentration Dosage

1.0%-10.0%

Storage

The product is liquid form, and can be stored for 18 months at 2~8°C in a sealed and dark environment. To prevent microbial contamination after opening, try to keep the cap tightly sealed, store carefully once opened, and use it up as soon as possible.

Mobel Biomaterials Technology Limited

Factory Address: 25 No. Bio-chemicals Technical Zoon, Yuhang District, Hangzhou city, Zhejiang P.R. China 310001

Homepage: www.mobelbiochem.com Email: infor@mobelhealth.com Tel: +86-571-85099223 Fax: +86-571-85099222